


Welcome to **navify**® Lab Operations

Simple. Seamless.

As healthcare becomes more complex, your lab works to make finding answers simple.

navify® Lab Operations is an advanced middleware solution that ensures seamless integration and operational efficiency so you can do more with less. Designed to work with the **cobas**® family of analyzers and automation systems, **navify**® Lab Operations delivers standardization and control to drive maximum performance in labs of all sizes.





Next-generation middleware that is anything but middle-of-the-road

navify® Lab Operations dynamically elevates the performance of your lab—regardless of size—by collating and analyzing data for actionable insights. This next-generation middleware enables patient data to be accessed securely anytime, anywhere for ultimate flexibility while allowing tighter integration among your automation, analytics, and IT systems. **navify**® Lab Operations dynamically adjusts to accommodate changing conditions so your lab can respond quickly to challenges throughout the day.

navify® Lab Operations was designed to simplify

- Its efficient workflow engine helps you **manage tasks and work processes** in your lab work area
- Dynamic routing capabilities automatically send your sample to the optimal route, **improving turnaround time** and **saving manual input**



Staffing issues may be unavoidable, but **navify**® Lab Operations helps bridge the gaps

- With its automated capabilities, **navify**® Lab Operations helps **maintain the turnaround time consistency and high-quality standards** your clinicians and patients need
- Easy-to-configure autoverification settings simplify lab processes to **expand the capabilities of your staff**—less complex rule writing and less additional manual labor for verifying results
- Sample tracking and archiving gives your lab **real-time information** about tube location, route, image, and detailed route history of both primary and aliquot specimens



Patients are counting on reliable results—**navify**® Lab Operations helps you locate them faster

- **Easily view the information you need on a single screen**, including test results, instrument alarms, tube image,* patient history, and more
- An intuitive interface for **minimal staff training** and **minimal staff downtime**
- Configure monitoring and sample validation screen to **focus on priority tasks**



When it comes to quality assurance, **navify**® Lab Operations reduces hands-on time to minimize errors

- Automatically stops transmission of patient results to LIS to **ensure quality standards when the analyzer is impacted by quality control issues**
- Integrated QC process that **saves resources and reduces downtime**
- Advanced 24/7 tracking of average of normals in real time so every sample is a quality check to ensure **confidence with instrument performance**



Consistency. Connectivity. Unmatched compatibility. **navify**® Lab Operations checks all the boxes

- Part of a comprehensive range of laboratory solutions designed to **seamlessly support your lab's growth**
- Simple and secure browser-based access for **connectivity across your network and compatibility with your existing IT infrastructure**
- Seamlessly integrate with automation and analytics—so your lab can **process more samples in less** time without sacrificing analytical quality


The quality you expect from Roche

For more than 125 years, Roche has built a reputation for excellence in instruments and healthcare expertise.

We're a healthcare company first, so when we develop our all-in-one solutions, **we start with your system's healthcare goals and build from the ground up.** Our attention to detail, accuracy, timeliness, and standardization ensure that solutions like **navify**® Lab Operations live up to our standards—and yours.



Leverage these capabilities to realize your lab's full potential



Workflow engine – reduce lab stress and improve turnaround time with automated decision-making



Autoverification – reduce complexity with easy-to-configure settings that simplify laboratory processes



Sample tracking and archiving – always know the tube location and route history, and easily locate samples in storage



Single-screen results viewing – get a comprehensive, integrated view of critical data using a streamlined interface



Average of normals – review patient data over time to respond to issues before they occur, and preserve result integrity



Integrated QC – monitor QC trends, compare your laboratories, and easily display affected/associated patient results after QC errors



Connectivity – use simple and secure browser-based access for connectivity and compatibility across your network



Dashboard and analytics – gain control of lab operations with customizable dashboards you can view at a glance

